

# GLOVE

## Permeation System with PTFE cells In-organic Chemicals

### APPLICABLE TO EN 16523-1

Gets your in-organic chemical permeation testing equipment from SATRA Technology, UK, the global known testing services in the world. The equipment is definitely the best solutions to meet your requirement. The system has been developed and used by SATRA to perform EN 16523-1 chemical permeation testing of gloves. KEMAS KURNIA is the exclusive distributor in providing after-sales service for the system in Malaysia.

DISCOUNT 15% ON  
SATRA  
PERMEATION  
SYSTEM!!  
Contact us ASAP

**In-depth knowledge from the Expert is our added value services to your of investment.**

When you equip your laboratory with SATRA's test equipment you will definitely obtain knowledge from the expertise of the field. With the assistance of SATRA's testing team, you will learn on how to specify a test procedure accordance to EN16523-1:2015. An on-site installation and training will boost your confidence in using the instrument along with the test process method. Assistance includes software analysis and interpretation of the testing result. **NEVER SETTLE FOR LESS** when buying an equipment, you need to **GO FOR MORE.**

**THINK SAFETY FIRST!**

SATRA never used glass cells in their chemicals permeation testing analysis. This is because glass is not a very safe material to use with chemicals. Glass is a fragile material and it has tendency to break, it has short life, since it was so susceptible to chemicals damaging the inner path exposed to acids and bases and of course we care about your safety first. Broken glass is a physical hazard. Using PTFE Cells in your testing will eliminate the potential risk of injury during testing. **CHANGE IS GOOD.**

Contact us via [info@kemaskurnia.com](mailto:info@kemaskurnia.com) or +6016-3280601 for more details the scope of the supply of the permeation system.